

ASSEMBLY INSTRUCTIONS

**MODEL : 0407A
CHEST**




**(BOX 1)
BOX OF CHEST**

NO	ITEM	QTY
A	CHEST	1 pc
B	KD LEG LEFT	2pcs
C	KD LEG RIGHT	2pcs

**(BOX 1)
HARDWARE LIST**

NO	ITEM	QTY
1	Csk m4 x 28mm Rbw	12 pcs

STEP 1 : Pull the bottom drawer and take out the KD leg carton (B & C). Assemble the KD leg (B & C) with CSK Screw(1) and securing with screw driver as shown.



COMPLETED ASSEMBLY

NOTE : MUST TIGHTEN SCREWS **PERIODICALLY** WITH USE SCREW WILL BECOME LOOSE CHECK TIGHTNESS OF ALL SCREW EVERY 6-8 WEEKS.